

## **Optimization of pulse labelling techniques of plant shoots with 14CO2 or 13CO2 to study the belowground allocation of recently assimilated C**

Christophe Robin

## **To cite this version:**

Christophe Robin. Optimization of pulse labelling techniques of plant shoots with 14CO2 or 13CO2 to study the belowground allocation of recently assimilated C. Rhizosphere2 International Conference, Aug 2007, Montpellier, France., Rhizosphere 2: book of Abstrats, Session 5 (P-960), pp.87, Cutting edge approaches and methodologies. hal-04474768

## **HAL Id: hal-04474768 <https://hal.univ-lorraine.fr/hal-04474768v1>**

Submitted on 23 Feb 2024

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

**Procedure**: Prior to the labelling, the belowground compartment is hermetically separated from the labelling atmosphere with liquid silicon set between shoots and roots (Photo 1). Plants are set into a labelling chamber (Photo 2) whose climatic parameters are regulated (temperature, relative air humidity and light).

The <sup>14</sup>CO<sub>2</sub> (or <sup>13</sup>CO<sub>2</sub>) is generated by the addition of NaH<sup>14</sup>CO<sub>3</sub> (or NaH<sup>13</sup>CO<sub>3</sub>) in 1M lactic acid (Photo 3). CO<sub>2</sub> and <sup>14</sup>C specific activity of the atmosphere in the labelling chamber are regulated with a β radiations scintillation detector (Photo 4) coupled with an infra-red gas analyser (IRGA), mass flow controller for <sup>12</sup>CO<sub>2</sub> and a peristaltic pump for

<sup>14</sup>CO<sub>2</sub> from NaH<sup>14</sup>CO<sub>3</sub>. For <sup>13</sup>C labelling, <sup>12</sup>CO<sub>2</sub> and <sup>13</sup>CO<sub>2</sub> can be regulated separately with a dedicated IRGA. The required <sup>14</sup>C/<sup>12</sup>C or <sup>13</sup>C/<sup>12</sup>C ratio is reached within 10 min. after the begining of the labelling and remains constant, regulated according to the shoot demand (Fig. 2). The belowground compartment is continuously flushed with decarbonated air at the input. CO<sub>2</sub> evolved from this compartment at the output is either measured with an IRGA or trapped for measurement of respired <sup>12</sup>C, <sup>14</sup>C or <sup>13</sup>C (Photo 2).

At the end of labelling, air is forced to pass through a soda lime trap in order to decrease the <sup>14</sup>C/<sup>12</sup>C or <sup>13</sup>C/<sup>12</sup>C ratio in the chamber, while <sup>12</sup>CO<sub>2</sub> is continuously injected in the chamber to compensate for trapping and shoot uptake. As a part of the <sup>14</sup>C (or <sup>13</sup>C) respired by the shoots during the chase period can be re-fixed during the chase period, the <sup>14</sup>C (or <sup>13</sup>C) is continuously fixed in a soda lime trap. Thus, the CO<sub>2</sub> concentration is steady during all the procedure (at 598 ± 7.5 vpm in the example shown Fig. 1).

Atmospheric CO<sub>2</sub> can be regulated from very low concentrations up to 1000 vpm. Pulse labelling or continuous labelling is feasible, either with a single C isotope or for dual <sup>13</sup>C and <sup>14</sup>C labelling. *See poster number by Koller et al. showing an example of <sup>13</sup>C pulse labelling.* 

Introduction of the Chapter 3.2 "Element cycling and organic matter turnover". In "Handbook of Methods in rhizosphere research". Action COST 631 « Understanding and Modelling the plant-soil interactions in the rhizosphere environment ». 2007. Jörg Luster and Roger Finlay (Eds). pp 62-68. ISBN 3-905621-35-5.

Allard V., Robin C., Soussana J.F., Liefering M., Newton P.C.D., 2006. Short and long-term effects of elevated CO<sub>2</sub> on *Lolium perenne* rhizodeposition and its consequences on soil organic matter turnover and plant N yiel *Biology and Biochemistry* 38: 1178-1187.

Bazot S., Ulff L., Blum H., Nguyen C., Robin C., 2006. Effects of elevated CO<sub>2</sub> concentration on rhizodeposition from *Lolium perenne* grown on soil exposed to 9 years of CO<sub>2</sub> enrichment. Soil Biology and Biochemistry 38 736.

Photo 4. Counting of  $14CO<sub>2</sub>$   $\beta$ radiations within a 15 ml plastic scintillant cell

Rhiz Sphere 2 Montpellier, 26-31 August 2007

**Figure 2**: Specific activity (<sup>14</sup>C/<sup>12</sup>C, in raw counts in the graph) of the atmosphere fixed at 50 KBq mg<sup>-1</sup> C during labelling. Arrows represent the beginning and the end of the labelling period.<sup>14</sup>C is detected with a scintillation detector Radiomatic A500 TR Packard Instruments (one acquisition every two min.; photo 4) and CO<sub>2</sub> is measured with an IRGA.

Photo 3. Generation of  ${}^{14}CO_2$  (or  ${}^{14}CO_2$ )



Objective : to describe a procedure of plant shoots labelling with <sup>14</sup>CO<sub>2</sub> or <sup>13</sup>CO<sub>2</sub>, with minor disturbances of the plant functioning. This is a prerequisite of studies aiming to investigate short-term assimilate partitioning in the plant and rhizodeposition. For such labelling procedure, it is particularly crucial to maintain steady the CO<sub>2</sub> concentration and the specific activity (or isotopic enrichment) of the atmosphere surrounding shoots.

## **Related references** :

Photo 1. Airtight separation of shoot and roots. The belowground compartment is flushed with decarbonated air at a flow 50-120 ml min-1



