



# Analysis of phase transformation induced microtextures in ODS steels: in-situ hightemperature EBSD experiments versus crystallographic reconstructions

Maissa Fekih, L Germain, Denis Sornin, J Guyon, Nathalie Gey

## ► To cite this version:

Maissa Fekih, L Germain, Denis Sornin, J Guyon, Nathalie Gey. Analysis of phase transformation induced microtextures in ODS steels: in-situ hightemperature EBSD experiments versus crystallographic reconstructions. ICOTOM 20, Jun 2024, Metz, France. hal-04787676

HAL Id: hal-04787676

<https://hal.univ-lorraine.fr/hal-04787676v1>

Submitted on 17 Nov 2024

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# **Analysis of phase transformation induced microtextures in ODS steels: in-situ high-temperature EBSD experiments versus crystallographic reconstructions**

M. Fekih<sup>1</sup>, L. Germain<sup>1,2</sup>, D. Sornin<sup>3</sup>, J. Guyon<sup>1,2</sup>, N. Gey<sup>1,2</sup>

<sup>1</sup>Université de Lorraine, CNRS, Arts et Métiers Paris Tech, LEM3, F-57000 Metz, France

<sup>2</sup>Laboratory of Excellence on Design of Alloy Metals for Low-mAss Structures (DAMAS), Université de Lorraine, France, <sup>3</sup>Université Paris-Saclay, CEA, Service de Recherches en Matériaux et procédés Avancés, 91191, Gif-sur-Yvette, France

Oxide dispersion strengthened (ODS) steels are designed for advanced nuclear reactors to enhance the high-temperature characteristics of body-centered stainless steel. To understand the microstructure evolution after austenitic heat treatments, Anthony Durand [1, 2] has utilized various analysis tools but they still insufficient to explain the microstructure inheritance and the associated phase transformation mechanisms.

Nowadays, the development of compact hot stages compatible with EBSD set-up as well as the emergence of high-speed EBSD camera open a promising way to track direct microstructure evolutions induced by phase transformations [4, 5].

In this work, a heating module for in-situ EBSD observations coupled with a CMOS EBSD camera, was used to observe directly phase transformation evolutions in 10Cr and 12Cr steel with and without ODS with different initial in 10Cr and 12Cr steel with and without ODS with different initial microstrucutres [1, 2]. Additionally, indirect parent-phase reconstruction has also been applied. The in-situ system was able to capture the direct transformation of ferrite/austenite during heating and the martensitic transformation during cooling of both steels with and without ODS and prove the presence of non-transformed ferrite at 1000°C in ODS Steel. The results allow to discuss the influence of the initial microstructure on the preferential nucleation site for austenite. **We aim to give a feedback on advantages and drawbacks of in-situ HT EBSD analysis applied to ODS steels and its complement to indirect crystallographic reconstruction.**

- [1] Durand, A. Thèse de doctorat, Université Paris-Saclay., 2021.
- [2] Durand, A. et al. *Materialia* 16, mai 2021: 101066. <https://doi.org/10.1016/j.mtla.2021.101066>.
- [3] Germain, L. et al. *Acta Materialia* 60, n° 11 (2012): 4551- 62. <https://doi.org/10.1016/j.actamat.2012.04.034>.
- [4] T. Martinez Ostormujof. Thèse de doctorat. Université de Lorraine, 2022. English. (NNT: 2022LORR0358). (tel-04278963)
- [5] Fekih, M. et al. «High temperature EBSD experiment applied to in-situ observation of phase transformation in steels» to be submitted